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FDA\EMA Immune Cell Therapy Regulations

Publishers	File name	Time of publication
FDA	Chemistry, Manufacturing, and Control (CMC) Information for Human Gene Therapy Investigational New Drug Applications (INDs)	2020
FDA	Testing of Retroviral Vector-Based Human Gene Therapy Products for Replication Competent Retrovirus During Product Manufacture and Patient Follow-up	2020
FDA	Long Term Follow-Up After Administration of Human Gene Therapy Products	2020
FDA	Human Gene Therapy for Retinal Disorders	2020
FDA	Human Gene Therapy for Rare Diseases	2020
FDA	Human Gene Therapy for Hemophilia	2020
FDA	Recommendations for Microbial Vectors used for Gene Therapy	2016
FDA	Design and Analysis of Shedding Studies for Virus or Bacteria-Based Gene Therapy and Oncolytic Products	2015
FDA	Considerations for the Design of Early-Phase Clinical Trials of Cellular and Gene Therapy Products	2015
FDA	Preclinical Assessment of Investigational Cellular and Gene Therapy Products	2013
FDA	Potency Tests for Cellular and Gene Therapy Products	2011
EMA	Guideline on quality, non-clinical and clinical requirements for investigational advanced therapy medicinal products in clinical trials	2019
EMA	Guidelines on Good Manufacturing Practice specific to Advanced Therapy Medicinal Products	2017

Scan the QR code for the FDA/EMA guidelines package

Immune Cell Therapy Development Process (take CAR-T cell drugs as an example)

	Early drug discovery	Manufacture/Quality control	Non-clinical studies	Clinical research
Research Content	<ul style="list-style-type: none"> Targets identification and selection. Preparation and screening of scFv or sdAb (single domain antibody). Design and optimization of CAR structure. 	<ul style="list-style-type: none"> Selection of raw materials Preparation of Gene vector. Preparation of CAR-T cell therapy products. 	<ul style="list-style-type: none"> Selection of animal models. Pharmacodynamic studies. Pharmacokinetics. Non-clinical safety. 	<ul style="list-style-type: none"> Exploratory Clinical Trials Safety and tolerability. Assessment of activity in vivo. Confirmatory clinical trials. Clinical efficacy and safety.
Solutions and Products	50 + CAR-T Target Products: Suitable for multiple scenarios such as immunology, single-domain antibody or scFv screening, CAR affinity detection, and species cross-validation.	Raw materials: Anti-CD3/CD28 antibody-coupled magnetic beads, cytokines, GENES™Nuclease, and other products. For quality control release testing: fluorescent-labeled proteins	Pharmacokinetic study: Highly sensitive and specific fluorescent-labeled CAR-T target proteins, anti-idiotypic antibodies. For immunogenicity evaluation: anti-idiotypic antibodies, such as rabbit Anti-FMC63 scFv polydonal antibody and murine Neutralizing Anti-FMC63 scFv monoclonal antibody; ACROBiosystems can provide customized anti-idiotypic antibody services to meet customers' demands.	

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


The purity of Human CD19 (20-291), Fc Tag (Cat. No. **CD9-H5251**) was more than 95% and the molecular weight of this protein is around 140-160 kDa verified by HPLC-MALS.

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
Human CD19 (20-291), Fc Tag (Cat. No. **CD9-H5251**) captured on CM5 chip via Anti-Human IgG Fc antibodies surface, can bind FMC63 MAb (Mouse IgG2a) with an affinity constant of 0.17 nM as determined in a SPR assay (Biacore T200) (Routinely tested).

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
Binding activity of FITC-Labeled Human CD19, His Tag from two different vendors were evaluated in the ELISA analysis against FMC63 Mab. The result showed that ACRO's FITC-Labeled Human CD19, His Tag has a much higher binding activity than that of the other vendor.

Binding activity of FITC-Labeled Human CD19, His Tag from two different vendors were evaluated in the flow cytometry analysis against anti-CD19-CAR-293 cells. The result showed that ACRO's FITC-Labeled Human CD19, His Tag has a much higher binding activity than that of the other vendor.

Binding activity of Human CD19, His Tag from two different vendors were evaluated by SPR assay against FMC63 mAb. The result showed that ACRO's Human CD19, His Tag can bind FMC63 mAb with an affinity constant of 2.95 nM which is much higher than that of the other vendor.

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

Immobilized Human HLA-A*02:01 NY-ESO-1 (SLLMWITQC) Tetramer Protein (Cat. No. **HL1-H52E8**) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Anti-NY-ESO-1 Antibody, Human IgG1 with a linear range of 0.2-8 ng/mL (QC tested).

Anti-NY-ESO-1 antibody captured on CM5 chip via Anti-human IgG Fc antibodies surface can bind Human HLA-A0201 NY-ESO-1 (SLLMWITQC) Tetramer Protein (Cat. No. **HL1-H52E8**) with an affinity constant of 1.35 nM as determined in an SPR assay (Biacore 8K) (QC tested).

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

GMP Human IL-15 (Cat. No. **GMP-L15H13**) stimulates the proliferation of CTLL-2 cells. The EC50 for this effect is 1.004 ng/mL, corresponding to a specific activity of $> 0.8 \times 10^7$ IU/mg, which is calibrated against human IL-15 WHO International Standard (NIBSC code: 95/554).

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

GMP Human IL-15 (Cat. No. **GMP-L15H13**) is stable in undiluted samples at 25°C for 72 hours and freeze-thaw 3 times without performance reduction.

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GMP Human IL-15 (Cat. No. **GMP-L15H13**) is stable in undiluted samples at 4 °C for 14 days without performance reduction.

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
Bioactivity of three different lots of GMP Human IL-15 (Cat. No. *GMP-L15H13*) verified by cell-based assay, and the result shows very high batch-to-batch consistency.

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
The purified human T cells were activated using Anti-CD3/CD28 Antibody-coupled Magnetic Beads (Cat. No. *MBS-C001*) at a ratio of 1:1 beads-to-cells for 24 hours along with RPMI1640 supplemented with 10% of FBS. The negative control experiment was performed by adding the Negative Control Beads coupled HSA. Cells were fluorescently stained using PE-labeled anti-human CD25 antibody and labeled FITC anti-human CD69 antibody and analyzed by flow cytometry.

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The purified human T cells were stimulated using Anti-CD3/CD28 Antibody-coupled Magnetic Beads (Cat. No. *MBS-C001*) at a ratio of 1:1 beads-to-cells. Cells were expanded in a T cell culture medium supplemented with 4ng/mL of rIL-2 Protein (Acrobiosystems, Cat. No. *IL2-H4113*). Activated Cells were expanded for up to 13 days (A) with high cell viability (B).

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


Measured by its ability to cleave a targeted DNA substrate. Cas9 (Cat. No. CA9-S5149) achieves >90% substrate cleavage.

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



2x10⁵ of 293 CAR cells transfected with anti-BCMA-scFv were stained with 100 μ L of 1 μ g/mL of FITC-Labeled Human BCMA, His Tag (Cat. No. **BCA-HF2H3**) and negative control protein respectively (Fig. C and B), and non-transfected 293 cells were used as a control (Fig. A). FITC signals were used to evaluate the binding activity (QC tested).

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
The purity of FITC-Labeled Human BCMA, His Tag (Cat. No. **BCA-HF2H3**) was more than 90%, and the MW of this protein is around 24-34 kDa as verified by HPLC-MALS.

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The binding activity of FITC-Labeled Human BCMA protein from AcroBiosystems and a competing vendor was evaluated by FACS analysis. The result showed that ACRO's Star Staining FITC-Labeled Human BCMA (Cat. No. **BCA-HF2H3**) protein has a much higher binding activity as compared to competitors.

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

Non-specific binding to non-transduced PBMCs between FITC-Labeled Human BCMA Protein of Acro and competitor. 5e5 of non-transduced PBMCs were stained with FITC-Labeled Human BCMA Protein and anti-CD3 antibody, washed and then analyzed with FACS. PE signal was used to evaluate the expression of CD3+ T cells in non-transduced PBMCs, and FITC signal was used to evaluate the non-specific binding activity to non-transduced PBMCs.

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
The binding activity of the Human BCMA before and after FITC labeling was evaluated in the above FACS analysis. The result shows that FITC-Labeled BCMA (Cat. No. *BCA-HF2H3*) and unconjugated Human BCMA have similar levels of binding activity.

Binding affinity of the Human BCMA before and after FITC labeling was evaluated in the above SPR analysis (Biacore T200). The result shows that FITC-Labeled (Cat. No. *BCA-HF2H3*) and unconjugated Human BCMA, His Tag have almost the same level of affinity.

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

Binding activity of two different sizes (left) and three different lots (right) of FITC-Labeled Human BCMA (Cat. No. *BCA-HF2H3*) against anti-BCMA CAR-293 cells was evaluated by flow cytometry. The result shows very high batch-to-batch consistency.

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Comparison of FITC-Labeled Human BCMA Protein accelerate samples

Comparison of FITC-Labeled Human BCMA Protein Freeze-thaw samples

FITC-Labeled Human BCMA (Cat. No. *BCA-HF2H3*) kept at different concentrations at 25°C for 48 hours and undergoing multiple freeze-thaw cycles shows no stability issue or performance reduction.

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


The FACS analysis evaluated the binding activity of FITC-Labeled Human CD19 protein from two different vendors. The result shows that ACRO's Star Staining FITC-Labeled Human CD19 (Cat. No. CD9-HF2H3) protein has a much higher binding activity than other competitors'.

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
5e5 PBMCs were stained with FITC-Labeled Human CD19 (20-291), His Tag (Cat. No. CD9-HF2H3), and anti-CD3 antibody, washed and then analyzed with FACS. PE signal was used to evaluate the expression of CD3+ T cells in PBMCs. FITC signal was used to evaluate the nonspecific binding activity to PBMCs (QC tested).

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

1e6 of the Anti-CD19 CAR-293 cells were stained with 100 μ L of 1:50 dilution (2 μ L stock solution in 100 μ L FACS buffer) of PE-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. No. FM3-HPY53) and negative control antibody respectively. PE signal was used to evaluate the binding activity (QC tested).

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
5e5 PBMCs were stained with PE-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. FM3-HPY53), and anti-CD3 antibody, washed followed by FACS analysis. FITC signal was used to evaluate the expression of CD3+ T cells in PBMCs, and PE signal was used to evaluate the nonspecific binding activity to PBMCs.

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The binding activity of two different batches of PE-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. No. FM3-HPY53) against anti-CD19 CAR-293 cells was evaluated by flow cytometry. The result shows very high batch-to-batch consistency.

PE-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. No. FM3-HPY53) were tested in different concentrations at 37°C for 72 hours with repeated freeze-thaw cycles and found to be stable without performance reduction.

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


The purity of Rabbit Anti-Mouse FMC63 scFv Polyclonal Antibody (Cat. No. FM3-S93) was more than 90%, and the molecular weight of this protein is around 130-145 kDa, as verified by HPLC-MALS.

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
Immobilized FMC63 scFv at 2 µg/mL (100 µL/well) can bind Rabbit Anti-Mouse FMC63 scFv Polyclonal Antibody (Cat. No. FM3-S93) with a linear range of 0.098-0.78 ng/mL (QC tested).

Immobilized anti-CD19 antibody (Clone: FMC63) at 1 µg/mL, add increasing concentrations of Rabbit Anti-Mouse FMC63 scFv Polyclonal Antibody (Cat. No. FM3-S93) and then add Biotinylated anti-CD19 antibody (Clone: FMC63) at 2 µg/mL. Detection was performed using HRP-conjugated streptavidin with sensitivity of 78 ng/mL (QC tested).

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

The purity of Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Carrier-free) (Cat. No. FM3-Y45A1) was more than 95%, and the molecular weight of this protein is around 140-160kDa as verified by SEC-MALS.

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
ELISA analysis shows that the binding of Human CD19, Fc Tag (Cat. No. CD9-H5251) to FMC63 scFv, His Tag was inhibited by increasing concentration of Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Clone Y45). The concentration of Human CD19, Fc Tag used is 5 µg/mL (100 µL/well). The IC50 is 0.084 µg/mL (Routinely tested).

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Immobilized FMC63 scFv, His Tag at 2 µg/mL (100 µL/well) can bind Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Clone Y45) with a linear range of 1-19 ng/mL. Anti-DNP antibody, mouse IgG1 (Cat. No. *DNP-M1*) was used as an isotype control (QC tested).

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Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. No. *FM3-Y45*) captured on CM5 chip via anti-mouse antibodies surface can bind FMC63 scFv with an affinity constant of 1.08 nM as determined in an SPR assay.

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Anti-Idiotypic Antibody Development Service

One-stop services from antigen preparation to test kit developments

mAb

CAR-T

PK/ADA assay

BsAb

ADC

ACRO proteins supplied for free

Neutralizing

Non-neutralizing

Drug target complex

Please Scan the QR code for more information

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分子	カタログ番号	製品説明
IL-15	GMP-L15H13	GMP Human IL-15
IL-7	GMP-L07H14	GMP Human IL-7
IL-21	GMP-L21H18	GMP Human IL-21
CD3	GMP-MC0323	GMP Monoclonal Anti-Human CD3 Antibody (OKT3)

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分子	カタログ番号	製品説明
BCMA	BCA-HF2H3	FITC-Labeled Human BCMA / TNFRSF17 Protein, His TagStar Staining
CD19	CD9-HF2H3	FITC-Labeled Human CD19 (20-291) Protein, His TagStar Staining
Mesothelin	MSN-HF2H3	FITC-Labeled Human Mesothelin / MSLN (296-580) Protein, His TagStar Staining
Siglec-2	SI2-HF2H4	FITC-Labeled Human Siglec-2 / CD22 Protein, His TagStar Staining
FMC63	FM3-FY57P1	FITC-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Y45)Star Staining
CD19	CD9-HP2H5	PE-Labeled Human CD19 (20-291) Protein, His TagStar Staining
BCMA	BCA-HP2H7	PE-Labeled Human BCMA / TNFRSF17 Protein, His TagStar Staining
Mesothelin	MSN-HP2H8	PE-Labeled Human Mesothelin / MSLN (296-580) Protein, His TagStar Staining
CD19	CD9-HA2H9	APC-Labeled Human CD19 (20-291) Protein, His TagStar Staining
BCMA	BCA-HA2H4	APC-Labeled Human BCMA / TNFRSF17 Protein, His TagStar Staining
Mesothelin	MSN-HA2H6	APC-Labeled Human Mesothelin / MSLN (296-580) Protein, His TagStar Staining

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>>> Alexa Fluor 647 Label

BCMA	CD19	Mesothelin	Glypican 3	CD147	CD300e
HER2	Siglec-2	Siglec-3	PD-1	OX40	MUC-1
TSLP R					

>>> Alexa Fluor 555 Label

BCMA	CD19	Mesothelin	Glypican 3	NKG2D	FAP
Nectin-4	CD37	protein L			

>>> Alexa Fluor 488 Label

BCMA	CD19	Mesothelin	Glypican 3	CD147	CD300e
HER2	Siglec-2	Siglec-3	PD-1	OX40	MUC-1
TSLP R					

>>> FITC Label

Glypican 3	HER2	CD7	Siglec-3	CD30	EGF R
EGFRVIII	GUCY2C	NKG2D	FAP	Nectin-4	CD37
protein L					

>>> PE Label

Glypican 3	NKG2D	FAP	Nectin-4	CD37	protein L
CD147	CD300e	HER2	Siglec-2	Siglec-3	PD-1
OX40	MUC-1	TSLP R			

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カタログ番号	動物種	製品説明	用途
CD9-BV4324b	Mouse	Biotinylated Anti-CD19 scFv, Fc, Avitag™ (FMC63) (MALS verified)	ADA assay development
FM3-Y45A1	Mouse	Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Y45) (Carrier-free) (recommended for ADA assay) DMF Filed	
FM3-S93	Rabbit	Rabbit Anti-Mouse FMC63 scFv Polyclonal Antibody (recommended for ADA assay) (MALS verified)	

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カタログ番号	製品説明	用途
MBS-C001	Anti-CD3/CD28 Antibody-coupled Magnetic Beads	Cell activation and expansion

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分子	カタログ番号	製品説明
Cas9	CA9-S5149	NLS-Cas9 Nuclease
Cas12a	CAA-L5149	NLS-Cas12a Nuclease
BenzNuclease	BEE-N3116	GENIUS™ Nuclease DMF Filed
BenzNuclease	NUE-S5119	GENIUS™ Nuclease, premium grade
BenzNuclease	CRS-A016	GENIUS™ Nuclease ELISA Kit

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*Where Proteins and Innovation
Advance Biomedicine*

Copyright Statement

“

This material is copyrighted by the Company. All rights in this material are reserved by the Company. Unless otherwise indicated in writing, all material in this material is copyrighted by the Company. No part of this material may be copied, photocopied or reproduced in any form or redistributed to any other person or used in any other manner which infringes the Company's copyright without the prior written authorisation of the Company.

”

Her2 BAFFR LAG-3
 Fc Receptor Siglec-10
 Biotinylated Protein
 PD-L1 VEGF165 CD3 epsilon
 CD20 PD-1 BCMA
 CD27 PVRIG
 CD47 PSMA
 CFGL1 TFPI
 Siglec-15 Integrin
 CD24 CD3E & CD3D CD20
 CD19 FcRn PCSK9
 IL-2 R alpha
 CAR-T Target Protein
 Glypican 3 Integrin MICA PD-1
 FcRn ADA Service CD30 MICA CD3E & CD3G
 EGF R B7-H3 BCMA CD30 MICA CD3E & CD3G
 Integrin TIGIT TGF-beta 1
 4-1BB Siglec-15
 Biotinylated Protein
 CD20 CD200 GTR Nectin-4
 VEGF165 CD73 FGL1
 PCSK9 IgG1 Fc CD40 BCMA Nectin-4
 SIRP alpha ADA Service Immune Checkpoint Protein PD-L1
 Nectin-4 Biotinylated Protein CD3E & CD3D PSMA
 SPR /BLI analytical service IL-2