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


MHC-peptide complex can stably bind to T cell receptor after tetramerization.

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MHC monomers have a low avidity and fast dissociation rates when binding to the TCR.

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





Biotinylated Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Complex Protein (Cat. No. **HL1-H82E6**) on SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.


The purity of Biotinylated Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Complex Protein (Cat. No. **HL1-H82E6**) is more than 90% and the molecular weight of this protein is around 45-65 kDa verified by SEC-MALS.

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5e5 of anti-NY-ESO-1 TCR-293 were stained with 100 μ L of 1:25 dilution (4 μ L stock solution in 100 μ L FACS buffer) of PE-Labeled Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Tetramer Protein (Cat. No. **HL1-HP2E5**) and negative control protein respectively. PE signal was used to evaluate the binding activity (QC tested).





5e5 of anti-NY-ESO-1 TCR-293 were stained with 100 μ L of 1:25 dilution (4 μ L stock solution in 100 μ L FACS buffer) of PE-Labeled Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Tetramer Protein (Cat. No. **HL1-HP2E6**) and negative control protein respectively. PE signal was used to evaluate the binding activity (QC tested).

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

Performance of PE-labeled Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Tetramer from ACROBiosystems (Cat. No. **HL1-HP2E5**) and other brand were evaluated by flow cytometry under the manufacturer's recommended conditions. ACROBiosystems' PE-Labeled NY-ESO-1 Tetramer Protein shows 3x higher MFI signal strength than competitor M, revealing much better TCR binding.

Performance of PE-labeled Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Tetramer from ACROBiosystems (Cat. No. **HL1-HP2E6**) and other brand were evaluated by flow cytometry under the manufacturer's recommended conditions. ACROBiosystems' PE-Labeled NY-ESO-1 Tetramer Protein shows 2x higher MFI signal strength than competitor M, revealing much better TCR binding.

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
Immobilized Biotinylated Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Complex Protein (Cat. No. **HL1-H82E6**) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. **STN-N5116**) precoated (0.5 µg/well) plate can bind Anti-NY-ESO-1 Antibody, Human IgG1 with a linear range of 0.1-8 ng/mL (QC tested).

Immobilized Biotinylated Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Complex Protein (Cat. No. **HL1-H82E6**) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. **STN-N5116**) precoated (0.5 µg/well) plate can bind Anti-HLA class I Antibody, Human IgG1 (W6/32) with a linear range of 0.1-1 ng/mL (Routinely tested).

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

Anti-HLA class I Antibody, Human IgG1 (W6/32) captured on Protein A Chip can bind Human HLA-A*11:01&B2M&KRASG12D (VGDAGVGK) Complex Protein (Cat. No. **HLD-H52H4**) with an affinity constant of 1.06 nM as determined in a SPR assay (Biacore 8K) (Routinely tested).

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
The complex protein was verified stable after incubation under 37 degrees for 24 hours, and no significant dissociation of HLA-A*2402 and B2M (Cat. No. HLM-H82W4) was observed from HPLC-SEC analysis.

The complex protein was verified stable after incubation under 37 degrees for 48 hours, and no significant dissociation of HLA-A*0201 and B2M (Cat. No. HLM-H82W3) was observed from HPLC-SEC analysis.

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
Biotinylated Human HLA-A*0201 & B2M & NY-ESO-1-Accelerated stability test at 37°C

Biotinylated Human HLA-A*02:01&B2M Complex Protein (Peptide free, MALS verified) (Cat. No. HLM-H82W3) is stable in undiluted samples at 37°C for 48 h and freeze-throw cycles are stable without performance reduction.

Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Tetramer Protein-Accelerated stability test at 37°C

Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Tetramer Protein (MALS verified) (Cat. No. HLI-H52E8) is stable in undiluted samples at 37°C for 48 h and freeze-throw cycles are stable without performance reduction.

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HLA	Allele/ Target Peptide	製品番号	製品説明
HLA-A Complex (peptide free)	A*0201 & B2M	HLM-H52H5	Unconjugated Human HLA-A*02:01&B2M Complex Protein (Peptide free, MALS verified)
		HLM-H82W3	Biotinylated Human HLA-A*02:01&B2M Complex Protein (Peptide free, MALS verified)
	A*0301 & B2M	HLM-H82W6	Biotinylated Human HLA-A*03:01&B2M Complex Protein (Peptide free, MALS verified)
	A*1101 & B2M	HLM-H52H7	Unconjugated Human HLA-A*11:01&B2M Complex Protein (Peptide free, MALS verified)
		HLM-H82W5	Biotinylated Human HLA-A*11:01&B2M Complex Protein (Peptide free, MALS verified)
	A*2402 & B2M	HLM-H82W4	Biotinylated Human HLA-A*24:02&B2M Complex Protein (Peptide free, MALS verified)
HLA-C Complex (peptide free)	C*0102 & B2M	HLM-H82E8	Biotinylated Human HLA-C*01:02&B2M Complex Protein (Peptide free)
	C*0303 & B2M	HLM-H82E9	Biotinylated Human HLA-C*03:03&B2M Complex Protein (Peptide free)
HLA-DQ Complex (peptide free)	A1*03:02 DQB1*03:03	HL3-H5285	Human HLA-DQA1*03:02&DQB1*03:03 Complex Protein (Peptide free)
HLA-A*0201 & B2M-Peptide Complex	NY-ESO-1 (SLLMWITQC)	HL1-HP2E5	PE-Labeled Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Tetramer Protein
		HL1-H82E6	Biotinylated Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Complex Protein (Monomer, MALS verified)
		HL1-H52E8	Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQC) Tetramer Protein (MALS verified)
	NY-ESO-1 (SLLMWITQV)	HL1-HP2E6	PE-Labeled Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQV) Tetramer Protein
		HL1-H82E4	Biotinylated Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQV) Complex Protein (Monomer, MALS verified)
		HL1-H52E7	Human HLA-A*02:01&B2M&NY-ESO-1 (SLLMWITQV) Tetramer Protein (MALS verified)
	WT-1 (RMFPNAPYL)	HLW-H82E5	Biotinylated Human HLA-A*02:01&B2M&WT-1 (RMFPNAPYL) Complex Protein (Monomer, MALS verified)
		HLW-H52E6	Human HLA-A*02:01&B2M&WT-1 (RMFPNAPYL) Tetramer Protein (MALS verified)
	Glycipan 3 (FVGEFFTDV)	HLG-H82E4	Biotinylated Human HLA-A*02:01&B2M&Glycipan 3 (FVGEFFTDV) Complex Protein (Monomer, MALS verified)
		HLG-H52E5	Human HLA-A*02:01&B2M&Glycipan 3 (FVGEFFTDV) Tetramer Protein (MALS verified)
	PRAME (SLLQHLIGL)	HLE-H82E3	Biotinylated Human HLA-A*02:01&B2M&PRAME (SLLQHLIGL) Complex Protein (Monomer, MALS verified)
		HLE-H52H5	Human HLA-A*02:01&B2M&PRAME (SLLQHLIGL) Tetramer Protein (MALS verified)
	MAGE-A1 (KVLEYVIKV)	HLM-H82E1	Biotinylated Human HLA-A*02:01&B2M&MAGE-A1 (KVLEYVIKV) Complex Protein (Monomer, MALS verified)

HLA	Allele/ Target Peptide	製品番号	製品説明
HLA-A*0201 & B2M-Peptide Complex	MAGE-A2 (YLQLVFGIEV)	HLM-H82E6	Biotinylated Human HLA-A*02:01&B2M&MAGE-A2 (YLQLVFGIEV) Complex Protein (Monomer, MALS verified)
	MAGE-A3 (KVAELVHFL)	HLM-H82E3	Biotinylated Human HLA-A*02:01&B2M&MAGE-A3 (KVAELVHFL) Complex Protein (Monomer, MALS verified)
	MAGE-A4 (GVYDGREHTV)	HLM-H82E5	Biotinylated Human HLA-A*02:01&B2M&MAGE-A4 (GVYDGREHTV) Complex Protein (Monomer, MALS verified)
		HLM-H52H6	Unconjugated Human HLA-A*02:01&B2M&MAGE-A4 (GVYDGREHTV) Complex Protein (Monomer, MALS verified)
	MAGE-A4 (KVLHVVRV)	HLP-H82E5	Biotinylated Human HLA-A*02:01&B2M&MAGE-A4 (KVLHVVRV) Complex Protein (Monomer, MALS verified)
		HLK-H52E5	Human HLA-A*02:01&B2M&MAGE-A4 (KVLHVVRV) Tetramer Protein (MALS verified)
	GP100 (YLEPGPVTA)	HLG-H82E6	Biotinylated Human HLA-A*02:01&B2M&GP100 (YLEPGPVTA) Complex Protein (Monomer, MALS verified)
	RHAMM-R3 (ILSLELMKL)	HLR-H82E6	Biotinylated Human HLA-A*02:01&B2M&RHAMM-R3 (ILSLELMKL) Complex Protein (Monomer, MALS verified)
	AFP (FMNKFIYEI)	HLA-H82E3	Biotinylated Human HLA-A*02:01&B2M&AFP (FMNKFIYEI) Complex Protein (Monomer, MALS verified)
	EW (QLQNPSYDK)	HLE-H82E4	Biotinylated Human HLA-A*02:01&B2M&EW (QLQNPSYDK) Complex Protein (Monomer)
	CMV pp65 (NLVPMVATV)	HLC-H82E4	Biotinylated Human HLA-A*02:01&B2M&CMV pp65 (NLVPMVATV) Complex Protein (Monomer, MALS verified)
	HPV16-E6 (TIHDIILECV)	HLH-HP2H7	PE-Labeled Human HLA-A*02:01&B2M&HPV16-E6 (TIHDIILECV) Tetramer Protein
	HPV16-E7 (YMLDLQPET)	HL7-HP2H7	PE-Labeled Human HLA-A*02:01&B2M&HPV16-E7 (YMLDLQPET) Tetramer Protein
	HBV (FLLTKILTI)	HLT-H82E7	Biotinylated Human HLA-A*02:01&B2M&HBV (FLLTKILTI) Complex Protein (Monomer, MALS verified)
	HBV (FLLTRILTI)	HLT-H82E9	Biotinylated Human HLA-A*02:01&B2M&HBV (FLLTRILTI) Complex Protein (Monomer, MALS verified)
	HIV Gag (SLYNTVATL)	HLH-H82E8	Biotinylated Human HLA-A*02:01&B2M&HIV Gag (SLYNTVATL) Complex Protein (Monomer, MALS verified)
HLA-A*0301 & B2M-Peptide Complex	KRAS (VVGAGGVGK)	HLS-H82E5	Biotinylated Human HLA-A*03:01&B2M&KRAS (VVGAGGVGK) Complex Protein (Monomer)
		HLS-H52H4	Human HLA-A*03:01&B2M&KRAS (VVGAGGVGK) Complex Protein (Monomer)
	KRAS (VVGACGVGK)	HLC-H82E6	Biotinylated Human HLA-A*03:01&B2M&KRASG12C (VVGACGVGK) Complex Protein (Monomer)
		HLC-H52H4	Human HLA-A*03:01&B2M&KRASG12C (VVGACGVGK) Complex Protein (Monomer)
		HLC-H52H7	Human HLA-A*03:01&B2M&KRASG12C (VVGACGVGK) Tetramer Protein
HLA-A*1101 & B2M-Peptide Complex	KRAS (VVGAGGVGK)	HLS-H82E3	Biotinylated Human HLA-A*11:01&B2M&KRAS (VVGAGGVGK) Complex Protein (Monomer, MALS verified)
		HLS-H52H5	Human HLA-A*11:01&B2M&KRAS (VVGAGGVGK) Complex Protein (Monomer, MALS verified)

HLA	Allele/ Target Peptide	製品番号	製品説明
HLA-A*1101 & B2M-Peptide Complex	KRAS (VVGADGVGK)	HLS-HA2H9	Alexa Fluor 647-Labeled Human HLA-A*11:01&B2M&KRAS (VVGAGGVGK) Complex ProteinStar Staining (Monomer)
		HLD-HA2H7	Alexa Fluor 488-Labeled Human HLA-A*11:01&B2M&KRAS (VVGAGGVGK) Complex ProteinStar Staining (Monomer)
		HLD-H82E5	Biotinylated Human HLA-A*11:01&B2M&KRASG12D (VVGADGVGK) Complex Protein (Monomer, MALS verified)
		HLD-H52H4	Human HLA-A*11:01&B2M&KRASG12D (VVGADGVGK) Complex Protein (Monomer, MALS verified)
		HLD-H52H8	Human HLA-A*11:01&B2M&KRASG12D (VVGADGVGK) Tetramer Protein
	KRAS (VVGAVGVGK)	HLV-H82E6	Biotinylated Human HLA-A*11:01&B2M&KRASG12V (VVGAVGVGK) Complex Protein (Monomer, MALS verified)
		HLV-H52H3	Human HLA-A*11:01&B2M&KRASG12V (VVGAVGVGK) Complex Protein (Monomer, MALS verified)
		HLV-H52H8	Human HLA-A*11:01&B2M&KRASG12V (VVGAVGVGK) Tetramer Protein
HLA-C*07:02:01 & B2M-Peptide Complex	RYR (RYRPGTVAL)	HLR-H52W2	Biotinylated Human HLA-C*07:02:01&B2M&RYR (RYRPGTVAL) Complex Protein (Monomer)

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*Where Proteins and Innovation
Advance Biomedicine*

Copyright Statement

“

This material is copyrighted by the Company. All rights in this material are reserved by the Company. Unless otherwise indicated in writing, all material in this material is copyrighted by the Company. No part of this material may be copied, photocopied or reproduced in any form or redistributed to any other person or used in any other manner which infringes the Company's copyright without the prior written authorisation of the Company.

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