
Full Length Active Gallery

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

Full-length Claudin18.2-VLP (Cat. No. CL2-H52P7) has been observed under an electron microscope to ensure that it is assembled correctly.

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
The mean peak Radius of Human CCR5 Full Length Protein-VLP (Cat. No. CC5-H52P3) is 130-170 nm with more than 95% intensity as determined by dynamic light scattering (DLS).


The mean peak Radius of Human CXCR4 Full Length Protein-VLP (Cat. No. CX4-H5219) is 110-170 nm with more than 95% intensity as determined by dynamic light scattering (DLS).

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
Immobilized Human CCR5 Full Length Protein-VLP (Cat. No. **CC5-H52P3**) at 2 µg/mL (100 µL/well) can bind Anti-Human CCR5 Antibody, Mouse IgG2a with a linear range of 2-78 ng/mL (QC tested).


Immobilized Human CXCR4 Full Length Protein-VLP (Cat. No. **CX4-H5219**) at 5 µg/mL (100 µL/well) can bind Human ulocuplumab, Human IgG4 with a linear range of 1-31 ng/mL (QC tested).

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
Human Claudin-18.2 Full Length Protein-VLP (Cat. No. **CL2-H52P7**) captured on CM5 Chip via Anti-Claudin-18.2 antibody can bind Anti-Claudin-18.2 antibody with an affinity constant of 0.374 nM as determined in a SPR assay (Biacore T200).

Fluorescent Human Claudin-18.2 Full Length Protein-VLP (Cat. No. **CL2-HF218**) captured on CM5 Chip via Anti-Claudin-18.2 antibody can bind Anti-Claudin-18.2 antibody with an affinity constant of 0.433 nM as determined in a SPR assay (Biacore T200).

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2e5 of Anti-Claudin-18.2 CAR-293 cells were stained with 100 µL of 3 µg/mL of Fluorescent Human Claudin-18.2 Full Length Protein-VLP (Cat. No. **CL2-HF218**) and negative control protein respectively, FITC signals was used to evaluate the binding activity.

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


Biotinylated Human CD20 Full Length, His,Avitag (HEK293) (Cat. No. **CD0-H82E5**) on SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.


Human CD133, His Tag (Cat. No. **CD0-H82E5**) on SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.

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
Immobilized Rituximab at 2 µg/mL (100 µL/well) can bind Human CD20 Full Length, His Tag, HEK293 (SPR verified) (Cat. No. **CD0-H52H3**) with a linear range of 0.4-3 ng/mL (in presence of DDM and CHS) (QC tested).


Immobilized Biotinylated Human GPRC5D, His,Avitag&Flag Tag (Cat. No. **GPD-H82D6**) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. **STN-N5116**) precoated (0.5 µg/well) plate can bind Anti-GPRC5D Antibody, Human IgG4 with a linear range of 0.001-0.25 µg/mL (QC tested).

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Biotinylated Human CD20, His,Avitag (HEK293) (Cat. No. **CD0-H82E5**) captured on Biotin CAP-Series S Sensor Chip can bind Rituximab with an affinity constant of 1.73 nM as determined in a SPR assay (in presence of DDM and CHS) (Biacore T200)

Anti-GPRC5D Antibody captured on CM5 chip via Anti-human IgG Fc antibodies surface can bind Human GPRC5D, Flag,His Tag (Cat. No. **GPD-H52D3**) with an affinity constant of 71 nM as determined in a SPR assay (in presence of DDM and CHS) (Biacore 8K)

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


Biotinylated Human CD20 Full Length, His,Avitag (Cat.No.CD0-H82E3) on SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.


Human GPRC5D, Flag,His Tag (Nanodisc)(Cat. No. GPD-H52D4) on SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.

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
Immobilized Rituximab at 5 µg/mL (100 µL/well) can bind Biotinylated Human CD20 Full Length, His,Avitag (Cat. No. CD0-H82E3) with a linear range of 0.008-0.125 µg/mL (QC tested).


Immobilized Biotinylated Human GPRC5D, His,Avitag&Flag Tag (Nanodisc) (Cat. No. GPD-H82D4) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Anti-GPRC5D Antibody, Mouse IgG1 with a linear range of 0.002-0.313 µg/mL (QC tested).

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
Anti-GPRC5D Antibody captured on CM5 chip via Anti-human IgG Fc antibodies surface can bind Human GPRC5D, Flag,His Tag (Nanodisc) (Cat. No. GPD-H52D4) with an affinity constant of 0.179 µM as determined in a SPR assay (Biacore 8K).


Loaded Rituximab on Protein A Biosensor can bind with Biotinylated Human CD20 Full Length, His,Avitag (Cat. No. CD0-H82E3) an affinity constant of 31.3 nM as determined in BLI assay (ForteBio Octet Red96e).

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2e5 of CD20-CAR-293 cells transfected with anti-CD20-scFv were stained with 100 µL of 3 µg/mL of Human CD20 / MS4A1 Full Length Protein, His Tag (Nanodisc) (HEK293) (Cat. No. CD0-H52H1) and negative control protein respectively, washed and then followed by PE anti-His antibody and analyzed with FACS (QC tested).

2e5 of CD20-CAR-293 cells transfected with anti-CD20-scFv were stained with 100 µL of 3 µg/mL of Biotinylated Human CD20 Full Length, His,Avitag (Cat. No. CD0-H82E3) and negative control protein respectively, washed and then followed by PE-SA and analyzed with FACS (QC tested).

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テクノロジープラットフォーム	VLP	界面活性剤ミセル	ナノディスク
Difficulty	+	++	+++
Cycle	8-10 weeks	10-12 weeks	>16 weeks
Productivity	+++ (mg)	++	+ (µg)
Application	Immunization ELISA/SPR/BLI Cell-based assay CAR detection	Immunization ELISA/SPR/BLI	Immunization ELISA/SPR/BLI Cell-based assay CAR detection
Advantages	Membrane protein with very high immunopotency. Closes to native conformation	High purity	Membrane protein with high purity, non-cytotoxicity and close to native conformation
Disadvantages	Difficult to quantify	Toxic, low coating rate	Higher cost and more difficult to assemble

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分子	製品説明	用途	カタログ番号
CD20	Human CD20 / MS4A1 Full Length Protein, His Tag (HEK293) (SPR verified)	免疫 ELISA SPR/BLI	CD0-H52H3
	Biotinylated Human CD20 / MS4A1 Full Length Protein, His, Avitag™ (HEK293)		CD0-H82E5
	Cynomolgus CD20 / MS4A1 Full Length Protein, His Tag		CD0-C52H8
	Canine CD20 / MS4A1 Full Length Protein, His Tag		CD0-C52H5
Claudin18.2	Biotinylated Human Claudin-18.2 Protein, His,Avitag™ (active membrane protein)		CL2-H82E3
	Human Claudin-18.2 Protein, His,Twin-Strep Tag (active membrane protein)		CL2-H5587
	Mouse Claudin-18.2 Protein, His,Twin-Strep Tag (active membrane protein)		CL2-M5585
	Cynomolgus Claudin-18.2 Protein, His,Twin-Strep Tag (active membrane protein)		CL2-C5586
Claudin-18.1	Human Claudin-18.1 Protein, His,Twin-Strep Tag (active membrane protein)		CL1-H5588
Claudin-9	Human Claudin-9 / CLDN9 Protein, His,Twin-Strep Tag (active membrane protein, SPR verified)		CL9-H5586
STEAP1	Human STEAP1 Full Length Protein, His Tag (HEK293)		ST1-H52H3
CD133	Human CD133 Protein, His Tag		CD3-H52H4
GPCR5D	Human GPCR5D Protein, Flag,His Tag		GPD-H52D3
	Biotinylated Human GPCR5D Protein, His,Avitag™&Flag Tag		GPD-H82D6
	Cynomolgus GPCR5D Protein, Flag,His Tag	GPD-C52D3	
CCR5	Human CCR5 Protein, Flag,His Tag	CC5-H52D1	
Buffer	200x DDM CHS buffer	膜タンパク質の溶解性と活性を維持する膜タンパク質または抗体を希釈するSPRまたはBLIアッセイでランニングバッファを準備する	DC-11

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