
GMP Grade 移行

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


Evaluation of purity by SDS-PAGE under reducing conditions results in a purity greater than 90% for GMP IL-15 (Cat. No. GMP-L15H13).

Evaluation of CTL-2 cell proliferation activity using GMP IL-15 (Cat. No. GMP-L15H13). Calculated EC50 is 1.004 ng/mL, corresponding to a specific activity greater than 0.8×10^7 IU/mg, calibrated against the relevant WHO International Standard. (NIBSC code: 9095/554).

Evaluation of bioactivity between GMP Human IL-15 (GMP-L15H13) lots. Cell-based assays were performed to evaluate bioactivity for each lot at varying concentrations, resulting in high batch-to-batch consistency.

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


Monoclonal anti-human CD3 antibody, Mouse IgG2a (Clone: OKT3) evaluation of purity by SDS-PAGE under reducing (R) and non-reducing (NR) conditions (insert). Purity was calculated to be greater than 95%. SEC-HPLC was also performed to confirm purity greater than 95%.

Evaluation of activation activity on purified human T cells. Anti-CD3/Cd28 antibody-coupled magnetic beads (Cat. No. MBS-C001) at a ratio of 1:1 (beads/cells) were cultured together for 24 hours with RPMI1640 supplemented with 10% of FBS. Negative control was performed by the addition of negative control beads coupled with HSA. Cells were fluorescently stained using PE-labelled anti-human CD25 and FITC-labelled anti-human CD69 antibodies and analyzed through flow cytometry.

Cell viability was evaluated by using Anti-CD3/CD28 magnetic beads to stimulate purified donor human T cells at a ratio of 1:1 (beads/cells). T cell culture medium supplemented with 4 ng/mL of IL-2 protein (Cat. No. IL2-H4113) was used for cell expansion up to 13 days with high levels of cell viability throughout the experiment.

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


Evaluation of purity by SDS-PAGE under reducing conditions results in a purity greater than 95% for Cas12a Nuclease (Cat. No. CAA-L5149).

Evaluation of Cas12a bioactivity by cleavage of target DNA substrate. Cas12a achieves a >90% substrate cleavage. (QC tested).

Evaluation of structural homogeneity and purity by SEC-MALS. Purity of Cas12a Nuclease was evaluated to be more than 90% with a molecular weight of around 135-165kDa.

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


Evaluation of purity by SDS-PAGE under reducing conditions results in a purity greater than 95% for GENIUS™ Nuclease (Cat. No. NUE-S5119).

Evaluation of specific activity for GENIUS™ Nuclease was found to be greater than 1.2×10^6 units/mg protein. One unit will completely digest sonicated salmon sperm DNA.

Evaluation of structural homogeneity and purity by SEC-HPLC. Purity of GENIUS™ Nuclease was evaluated to be more than 99%.

BIOSYSTEMS
Acro

Key raw materials for cell and gene therapy

Cytokines

Antibodies

Activation Magnetic Beads

CAR/TCR-T Targets

GENIUS™ Nuclease

CRISPR-Cas Nuclease

<< GMP and Premium grade available
<< Support CMC manufacture and clinical research

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原材料の特性	内容
Material Identification & Microbial or Viral Contamination-free	<ul style="list-style-type: none">■ Molecular composition or formulation■ Health statement, TSE/BSE-free, pathogen testing■ Abnormal / Acute toxicity tests■ Required viral testing and screening documentation for materials utilizing human-derived materials.
Purity Evaluation	<ul style="list-style-type: none">■ Assays to detect and identify residuals■ Identification and documentation of impurities■ Documentation of purity and purity evaluation
Batch-to-batch Consistency	<ul style="list-style-type: none">■ Certificate of Analysis (CoA) provided by raw material supplier determines the lot-to-lot consistency.■ GMP-manufactured materials inherently demonstrate consistency.
Storage & Stability	<ul style="list-style-type: none">■ Recommended storage conditions (e.g. temperature, light, humidity) demonstrate consistent performance.■ Stability test results that outline product shelf-life are different from raw material products.

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	プレミアム Grade 製品	GMP Grade 製品
Application	Research and Development; Preclinical research to transition into clinical phases.	Designed to meet clinical phase requirements and bolster your IND application to various regulatory bodies.
Quality System	ISO 9001 and 13485 Quality Management System	ISO 9001 and 13485 Quality Management System (R&D Stage) including GMP Quality System (Production Stage)
Production	ISO-certified facilities	GMP-grade certified facilities
	Animal-origin free materials & BSE/TSE-free	Animal-origin free materials & BSE/TSE-free
	Transient or stable cell lines	Stable cell lines (Full inspection by USP / ChP)
	Pharmaceutical-grade materials	Pharmaceutical-grade materials
	Strict secondary sterile filtration	Strict secondary sterile filtration
	Class C+A cleanroom with manual fill-finish	Class B+A cleanrooms with automated fill-finish
	No additional viral clearance steps	2 additional viral removal and inactivation steps (nanofiltration + low pH)
Quality Control	Sterility / Mycoplasma testing	Sterility / Mycoplasma testing
	Endotoxin control and detection	Endotoxin control and detection
	Validated key production equipment and analytical instrumentation	Validated equipment, analytical instrumentation, experimental methods (audit trail available)
	Residual DNA/HCP testing	Residual DNA/HCP testing
	No additional quality control tests	Complete adventitious agent testing (virus testing and animal in vivo safety experiments)
Documentation	Common regulatory support	Comprehensive Regulatory Support Files
	DMF Files (Few products)	DMF Files (ALL products)

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